








SUPPLIES & ACCESSORIES CARDIAC OUTPUT



For use with Spacelabs Command Modules (90496 & 91496) configured with the cardiac output option C.



Order Code	Description	Length/Size	Qty
Cardiac Output Cables & Probes			
Cables			
	700-0027-00 TruLink Cardiac Output Adapter Cable	30 cm/12 in.	1 Each
	306655-003 Cardiac Output cable compatible with Baxter Edwards/Abbott injectate temperature probes and most models of thermodilution catheters		1 Each
	306655-004 Cardiac Output cable compatible with Becton Dickinson/Datex-Ohmeda injectate temperature probe		1 Each
Injectate Probes			
	690-0023-00 In-line injectate temperature probe Compatible with Datex-Ohmeda SP4045		1 Each
	690-0024-00 Reference solution in-line injectate temperature probe. Compatible with Datex-Ohmeda SP5030		1 Each
	690-0033-00 In-line injectate temperature probe. Compatible with Datex-Ohmeda SP4042. Not compatible with Abbott or Baxter in-line injectate probes or delivery systems.		1 Each
	690-0043-01 In-line injectate sensor housing, single use, sterile. Compatible with Datex-Ohmeda SP5045		25/Box
Closed Loop Injectate Kits			
	690-0031-00 CritiKit SP4572 CO in-line room temperature injectate kit, closed-loop delivery system, single use, sterile		5/Box
	690-0032-00 CritiKit SP4500 CO in-line ice temperature injectate kit, closed-loop delivery system, single use, sterile		4/Box