

Trolleys



O10-2131-00



MT-801



MT-201

SL6A Options

- 044-0008-00 FDA-XML Export
- 044-0009-00 DICOM Export

SL12A Options

- 044-0010-00 FDA-XML Export
- 044-0011-00 DICOM Export
- 044-0012-00 VCG



Resting ECG Recorders
6 & 12 Channel



CardioExpress SL6A



- Portable, lightweight design
- Convenient for mobile use and offices with limited space

- Display
- 5.6" Color TFT LCD Touch Screen

- Easy data input and operation
- Alphanumeric keyboard and one-touch operation

- Power supply
- Built-in rechargeable battery, AC/DC power supply

- Pacemaker detection
- Pacemaker detection that meets the requirements of ANSI/AAMI EC11

- Interpretation
- Automatic measurement and interpretation for adults and pediatrics

- R-R analysis
- The R-R analysis mode calculates related indicators of HRV (Heart Rate Variability).

- Recording
- Internal thermal printer and external printer (contact Spacelabs for a list of compatible printers)

- Memory and data transmission
- Supports external archiving to the Spacelabs Healthcare Sentinel System

- Wi-Fi
- Built in Wi-Fi

- Data Storage
- Internal: 500 10-second 12-channel recordings
 - External 4Gb: 20000 for DAT format and <20000 for SCP format

- Orders
- Supports orders for patient tests to be taken

Options

- Sentinel Cardiology Information Management System
- 1D & 2D Barcode readers
- Roll stand and carrying bag

CardioExpress Connectivity



Today's hospitals require a cardiology data management system that will facilitate analysis and reporting while also improving efficiency and eliminating errors. Our flexible and scalable **CARDIOLOGY INFORMATION MANAGEMENT SYSTEM (CIMS)**, Sentinel offers that and more – making the most of your product and IT investments while helping you to improve patient care.

Sentinel fully integrates reports from our entire range of non-invasive diagnostic cardiology procedures: Holter monitoring, ambulatory blood pressure monitoring (ABPM), 12-lead ECG, and ECG stress testing. Sentinel also provides an interface to a number of other manufacturers' ECG devices, enabling Sentinel to provide a cost effective solution to integrate data from multiple devices.

With Sentinel, you have immediate access to detailed information for a given patient, including historic ECG data that can also be integrated with your EMR/HIS system. Patient details and recording data is organized and stored in a central, networked database which reduces search times and improves efficiency.

We designed Sentinel to be highly configurable, enabling you to work the way you want to. When you control the workflow, you have greater flexibility to coordinate your procedures, further improving your day-to-day efficiency.

CardioExpress SL12A

CardioExpress SL12A is a twelve-channel ECG designed for high volume use in busy departments. It features advanced design and ergonomics to enable the operator to quickly and easily control the unit to meet individual requirements.



Features

- Compact design
- 12.1 inch foldable high resolution color touchscreen

- Easy data input and operation
- Alphanumeric keyboard

- Power supply
- Built-in rechargeable Li/ion battery

- Pacemaker detection
- Pacemaker detection that meets the requirements of ANSI/AAMI EC11

- Interpretation
- Automatic measurement and interpretation tested with authoritative CSE database
 - Complete digital filters, minimising baseline drift, AC and EMG interference
 - Real-time waveform freeze
 - 120 second ECG waveform review
 - Waveforms and diagnosis results preview before printing
 - Automatic arrhythmia detection and recording

- R-R analysis
- The R-R analysis mode calculates related indicators of HRV (Heart Rate Variability).

- VCG (Optional)
- VCG (Vectorcardiogram)* to show phase shifts of heart activity.

- Memory and data transmission
- Supports external archiving to the Spacelabs Healthcare Sentinel System

- Data Storage
- Internal: 500 10-second 12-channel recordings
 - External 4Gb: 20000 for DAT format and <20000 for SCP format

Options

- Internal storage of up to 500 10sec Resting ECGs
- Sentinel Cardiology Information Management System
- Barcode Reader
- Roll stand and carrying bag

* **Note:** VCG reports cannot be saved to the electrocardiograph or transmitted to Sentinel.